

# Sample Preparation Tutorial

---

Srinivas Chakravarthy

BioCAT / Sector 18ID, APS

Illinois Institute of Technology

# SAXS – The Basic Idea

## A PRELIMINARY CHARACTERISATION


### PURITY



SDS-PAGE and GEL FILTRATION

High molecular weight species must be removed


### CONCENTRATION



MEASURE ABSORBANCE

Also use A280:A260 to detect nucleic acids

### MONODISPERSITY



DYNAMIC LIGHT SCATTERING

Sample must be aggregate-free


### MATCHED SOLVENT



DIALYZE SAMPLE

Matched solvent required for accurate solvent subtraction

## B SMALL-ANGLE SCATTERING




DATA REDUCTION

SOLVENT SUBTRATCTION


$I(q)$

$$q = \frac{4\pi \sin \theta}{\lambda}$$



FOURIER TRANSFORM

$P(r)$




$$R_g^2 = \frac{\int P(r)r^2 dr}{2 \int P(r) dr}, \quad I(0) = 4\pi \int_0^{D_{max}} P(r) dr.$$

# SAXS – The Basic Idea

## C DATA VALIDATION

### INITIAL INSPECTION



#### GUINIER PLOT

Linearity is necessary but not sufficient for further analysis

### RADIATION DAMAGE



#### MEASURE TIMECOURSE

$I(0)$  and  $R_g$  should show no time-dependence

### CONCENTRATION EFFECTS



#### MEASURE CONCENTRATION SERIES

$I(0)/C$  and  $R_g$  should remain constant. Increase with  $C$  indicates aggregation, decrease indicates interparticle interference

### STANDARDS




#### MEASURE WATER AND/OR SECONDARY PROTEIN STANDARD

Normalise scattering to standard. Calculated molecular weights will be accurate for high quality samples

$$I(0) = N(\Delta\rho V)^2 = \frac{C\Delta\rho^2 v^2 MW}{N_A}$$

## D MODELING

### HIGH-RESOLUTION VALIDATION



#### CRYSOLO/CRYSON

Compare atomic models to solution scattering


### MULTI-DOMAIN MODELLING



#### DAMMIN/DAMMIF/SASREF/BUNCH

Ab initio and/or rigid-body refinement against scattering data

### SUBUNIT MODELING



#### MULTICH/MONSA/SASREF7

Use neutron contrast variation data to extract component scattering functions, and/or model subunits

### FLEXIBILITY



#### ECM

Model an ensemble of rigid structures against the scattering data

# Most Encountered Sample Quality Issues

---

- Concentration effects (interparticle interaction or repulsion)
- Radiation Damage (Aggregation or Degradation)
- Incorrect Buffer Subtraction
- Polydisperse Sample

# Preparation, Preparation, Preparation

---


- Arguably the most important part of a SAXS experiment is sample preparation.
- Optimization of the purification protocol.
- Biochemical / Biophysical characterization.
- Informed decisions towards the most suitable constructs (truncations / fusions).
- Talk extensively with the beam line staff before the trip about the experiment - eliminate guess-work from the experiment setup.

# Diagnostics


---

- SDS-PAGE gels and native gels – Be wary though, no obvious contaminants does not always translate to pristine SAXS data.
- Symmetric single peaks obtained by size exclusion chromatography – But there may be shelf life issues and storage and shipping related complications.
- No visible particulates and precipitate and no significant loss of sample during concentration.
- DLS – often not very reliable.
- SEC-MALS and ultracentrifugation.

# Size-Exclusion Chromatography In-Line with Multiple Angle Light Scattering





**SEC-MALS Instrumentation Schematic**  
Courtesy: Kushol Gupta (Univ. of Pennsylvania)



SEC-MALS with a nucleosome core particle  
Here demonstrates the need for assessing the quality of sample using a variety of techniques.



# SEC-MALS and Analytical Ultracentrifugation as Complementary Techniques








Molecule or complex	SEC-MALS	
	Molecular weight (kDa)	
	Observed (%)	Calculated
BSA	66.47 ± 2	67.0
MeCP2 full length	56.00 ± 4	52.44
A-Nuc147	198.0 ± 2	205.0
A-Nuc147-MeCP2 full length	251.4 ± 2	257.9 (1:1)
W-Nuc165	202.8 ± 1	208.95
W-Nuc165-MeCP2 full length	253.9 ± 0.6	261.4 (1:1)





# Atomic Coordinates Compared to SAXS Data



# Optimal Sample Concentration

$$I(q) \propto \frac{N}{V} V_{particle}^2 (\rho_1 - \rho_2)^2 FF(q) S(q)$$

“solution part”  
“protein part”

- Scattering intensity proportional to the molecular mass and the concentration.
- Dilution series, extrapolation to zero concentration, important to exclude concentration effects.
- Useful rule of thumb to estimate optimal concentration to start with – Molecular mass (Kda) \* Concentration (mg/ml) ~ 100
- Generally 0.25-10 mg/ml (higher for smaller proteins and higher q-ranges, most beamlines use 30 – 100 ul sample volume.
- For LC-SAXS using a superdex-200 /75 (or a superose-3/6/12) 10/300 column (column volume = 24ml), we usually use ~ 200ul of 3-10mg/ml sample).
- Nucleic acids scatter more than proteins, so much better signal at lower concentrations.

# Sample Concentration Calculation

---

- Theoretical Extinction Coefficients can have upto 5-10% discrepancy.
- Experimental determination using denatured material (6M Guanidium HCl) - see ***Calculation of protein extinction coefficients from amino acid sequence data*** Stanley C. Gill, a and Peter H. von Hippel *Analytical Biochemistry* Volume 182, Issue 2, 1 November 1989, Pages 319-326.
- Bradford, Lowry, Biuret, and other colorimetric assays (use multiple concentrations and work with concentration-range relevant to the SAXS experiment).
- DNaseI digestion of polynucleotides.
- SAXS (IO) and SEC-MALS (RI).


# Considerations for Buffer Composition

---

- Bring abundant amounts of MATCHED buffer (10X stocks when possible recommended for convenience of shipping and the opportunity to trouble shoot on site.
- Variety of buffer systems are suitable (Tris, HEPES, PBS etc...)
- Moderate salt concentrations – more than 1M NaCl not recommended but in extreme cases, can be done.
- Glycerol helps with prevention of radiation damage but higher concentrations than ~ 5% make buffer subtraction difficult and also decrease contrast.
- Avoid detergents. Try using concentrations well below CMC if you absolutely need them.
- Free radical scavengers such as DTT (1mM) or other reducing agents (BME, TCEP etc...), highly recommended (use fresh stocks).
- Retrieve used sample if possible for analysis in order to make better informed decisions during subsequent experiments.

# Measures for Optimal Buffer Matching

---



Dialysis



SEC Buffer





Buffer Exchange

# Inter-particle Interactions and Preventive Measures


---

- Solvent exposed charged residues – electrostatic repulsion or attraction – can be neutralized by increasing salt in buffer.
- Hydrophobic interactions – can be ameliorated by moderate amounts (sub-CMC) of detergents.
- Disulphide bonds – can be reduced by reducing agents such as DTT and TCEP.

# Manifestation of Inter-particle Interactions in SAXS



No Inter-Particle  
Interference




Aggregated

Inter-Particle  
Repulsion



# Post-Aggregation First-Aid



Centrifugation



If the shelf life is really short, on-site repurification or LC-SAXS (if available).





Spin Filters

Syringe Filters



Ultracentrifugation (helps in some cases but may be excessive in most, especially large proteins and complexes)

# Radiation Damage






## Remedies:

1. Modify buffer composition – a. Add free radical scavengers such as DTT, and TCEP. b. Low percentages (1-5%) of glycerol.
2. Alter Data collection strategy – continuous uni-directional sample flow during exposure (need a lot of sample), reduce exposure time, temperature, and beam attenuation.

# Low-UV signal is not necessarily low scatter.

---







LC-SAXS Configuration  
Schematic

Mathew, E, et al., [J Synchrotron Radiat.](#)  
2004 Jul 1;11(Pt 4):314-8. Epub 2004 Jun 23.

Integrated Intensity of Individual  
Exposures Plotted in Sequence

# SEC-SAXS: Truly Efficient Way to Collect Bio-SAXS Data



# Easy Access to SEC-SAXS Logistics

Commercially Available Size-Exclusion Columns Suitable for SEC-SAXS				
Column Name	Exclusion Limit (M <sub>r</sub> )	Recommended Flow Rate (ml/min)	Theoretical Plates*	Column Volume (ml)
Superdex-200 Increase 10/300	~ 1.3 X 10 <sup>6</sup>	0.75	> 48,000 m <sup>-1</sup>	24
Superdex-200 Increase 5/150	~ 1.3 X 10 <sup>6</sup>	0.45	> 42,000 m <sup>-1</sup>	3
Superdex-75 5/150	~ 1 X 10 <sup>5</sup>	0.15	> 25,000 m <sup>-1</sup>	3
Superdex-75 10/300	~ 1 X 10 <sup>5</sup>	0.50	> 25,000 m <sup>-1</sup>	24
Superose-6 10/300	~ 4 X 10 <sup>7</sup>	0.50	> 30,000 m <sup>-1</sup>	24
Superose-6 5/150	~ 4 X 10 <sup>7</sup>	0.15	> 30,000 m <sup>-1</sup>	3
Superose-12 10/300	~ 2 X 10 <sup>6</sup>	0.50	> 40,000 m <sup>-1</sup>	24
KW 802.5	~ 1.5 X 10 <sup>5</sup>	1.00	> 21,000 per column	~ 15
KW 803	~ 7 X 10 <sup>5</sup>	1.00	> 21,000 per column	~ 15
KW 804	~ 1 X 10 <sup>6</sup>	1.00	> 16,000 per column	~ 15
KW 402.5-4F	~ 1.5 X 10 <sup>5</sup>	0.33	> 35,000 per column	~ 5
KW 403-4F	~ 6 X 10 <sup>5</sup>	0.33	> 35,000 per column	~ 5
KW 404-4F	~ 1 X 10 <sup>6</sup>	0.33	> 25,000 per column	~ 5
KW 405-4F	~ 2 X 10 <sup>7</sup>	0.33	> 25,000 per column	~ 5


The number of **theoretical plates** can be calculated from a chromatographic peak after elution.  $N = 5.55 tR^2/w_{21/2}$  .

tR = The time between sample injection and the peak reaching a detector at the end of the column is termed the *retention time*.


w<sub>1/2</sub> is the peak width at half-height.






# LC-SAXS Setup at BioCAT (APS)



# No SEC




# SEC-SAXS





# PSTK, SepSecS and tRNA<sup>Sec</sup> form a stable complex



Slide Courtesy: Simonović M,  
University of Illinois at Chicago

# SEC Instrumentation Upgrades to Facilitate Automation & Combined Techniques



- Acquired an AKTA-pure – modular in design. Therefore amenable to adding elements in the future that will enable automatic sample loading (column valve, loop valve, etc...).
- We will interface with auto-sampler, MALS system from Wyatt Scientific and refractometer etc..
- Already has multiple wavelength UV detector that is capable of measuring absorbance at as many as 3 wavelengths simultaneously.
- Also made available smaller columns which can be used where appropriate to increase efficiency.

# LC-SAXS : Advantages

---

- If right column is chosen, proteins / complexes of interest can be separated from contaminants, degradation products, and/or aggregates immediately before exposure to X-rays.
- The elution peak, if symmetric, and if  $R_g$  remains unchanged, represents a dilution series and therefore extrapolation to zero concentration is achieved easily.
- Radiation damage is prevented due to unidirectional flow of sample therefore preventing repeated exposure.
- Buffer blank is obtained quite easily by averaging the regions before and after the peak.